



Trimount™

Product Description

Trimount™ is a 50gsm white mounting neutral pH tissue coated on both sides with a heat activated solvent acrylic adhesive, which is tack-free at normal ambient temperature and therefore does not require a release liner
Available in various sizes

Suggested Applications

Trimount™ is a dry mounting tissue for use in all makes of all dry mounting presses. Trimount™ is recommended for bonding RC photographic paper and printed materials to smooth surfaces.

Product Structure

Adhesive	Solvent Acrylic (Heat Activated)
Carrier:	PAPER / OFF WHITE
Adhesive	Solvent Acrylic (Heat Activated)

General Properties

Film Thickness	50gsm
Adhesive Layer	12gsm
Ratio Film/Adhesive	2/1

Process Settings

Roller Laminator NOT RECOMMENDED	Press Temperature range: 85°C to 100°C (185°F to 212°F) Dwell time: 30 to 90 minutes
--------------------------------------------	---------------------------------------------------------------------------------------------------

Specifications subject to change without notice.

Published information concerning Drytac products is based upon research which the company believes to be reliable and accurate. Such information is intended only as a guide and is given without guarantees. Purchasers should independently determine, prior to use, the suitability of each material for their specific purpose. Drytac does not accept liability for risks associated with the performance of the materials in application. All Drytac products are subjected to carefully quality control throughout the manufacturing process and are sold subject to Drytac standard conditions of sale

Drytac Europe Ltd (UK) – Drytac House, Filwood Road, Fishponds, Bristol, BS16 3RY, UK. Tel: +44 (0) 845-070-0660
Drytac Corp. (USA) – 5383 Glen Alden Drive, Richmond, Virginia 23231, USA. Telephone: 1-800-603-6610
Drytac Canada Inc. (Canada) – 220 Caldari Road, Concord, Ontario, L4K 3X5, Canada. Telephone: 1-800-353-288

